Pump In Style® Advanced

EN Instructions for use | ES Instrucciones de uso | FR Mode d'emploi

#101036449

#101036445

#101036453

PLEASE READ AND SAVE THESE INSTRUCTIONS. FOR MOST CURRENT BREAST PUMP INSTRUCTIONS VISIT WWW.MEDELABREASTFEEDINGUS.COM

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Important Safeguards

When using electrical products, especially when children are present, basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

DANGER: To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place or drop into water or other liquid.

⚠ WARNING: To avoid fire, electrocution, or serious burns:

- Do not leave product unattended when plugged into an electrical outlet.
- Always unplug electrical product immediately after use.
- Close supervision is necessary when this product is used near children or persons with disabilities to prevent strangulation by the tubing or power adaptor cord.
- Close supervision is necessary when this product is used near children or persons with disabilities to prevent choking on small parts.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep cord away from heated surface.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors, or operate where aerosol (spray) products are being used or where oxygen is being administered.
- Always inspect power adaptor and battery pack wires prior to use for damage or exposed wire. If damage is found, immediately discontinue use of power adaptor or battery pack and call Medela Customer Service at 1-800-435-8316.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- The breast pump and accessories are not heat-resistant: keep away from heated surfaces or open flames.
- Do not use near flammable materials.
- Do not use an electrical outlet device that has been exposed to water or other liquids including:
 - Do not use while bathing or showering.
 - Never place or drop into water or other liquids.
 - Do not run water over breast pump.
 - If a device has been exposed to water or other liquids, do not touch, unplug the device from electrical outlet, turn off and contact manufacturer.

The defined signal words identify all instructions that are important to safety. Failure to observe these instructions can lead to injury or damage to the breast pump or yourself. When used in conjunction with the following words, the defined signal words stand for:

WARNING Can lead to serious injury or death.

CAUTION Can lead to minor injury.

NOTE Can lead to material damage.

INFO Useful or important information that is not related to safety.

WARNING: To avoid health risk and reduce the risk of injury:

- This product is intended for use by a single user only and should not be shared between users.
- Do not operate a moving vehicle while pumping.
- Inspect all appropriate pump components before each use.
- Do NOT continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- If tubing becomes moldy, discontinue use and replace tubing.
- If breast milk backup occurs, you must clean tubing, faceplate and diaphragm before your next pumping session.
- Use the product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never use while sleeping or when drowsy.
- This device cannot be serviced or repaired. Do not repair yourself. Do not modify the device.
- Never use a damaged device. Replace damaged or worn parts.
- Before each use visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- Use only the power adaptor that comes with the breast pump.
- Pumping can induce labor. Do not pump until after giving birth. If you become pregnant while breastfeeding or breast pumping, consult with a licensed healthcare professional before continuing.
- If irritation or discomfort occurs, discontinue use and see a doctor.
- If infected with Hepatitis B, Hepatitis C, or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- Do not microwave or boil breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)
- Clean and sanitize all parts that come into contact with your breast and breast milk prior to first use.
- Wash all parts that come into contact with your breast and breast milk after every use.
- Only use Medela recommended parts with your breast pump.
- Close supervision is necessary when the breast pump or accessories are used in the vicinity of children.

⚠ CAUTION: Can lead to minor injury:

- Do NOT wrap cord around adaptor body.
- Use only the power adaptor that comes with the Pump In Style® Advanced Breast Pump.
- Use only the Pump In Style Advanced battery pack that comes with the Pump In Style Breast Pump.
- Pump only with the breast pump in an upright position.
- Do NOT shorten tubing.
- Make sure the voltage of the power adaptor is compatible with the power source.
- Plug the power adaptor into the breast pump first and then into the wall socket.
- Do NOT use antibacterial or abrasive cleaners/detergents when cleaning breast pump or breast pump parts.
- Never put breast pump motor in water or a sterilizer, as you can cause permanent damage to the breast pump.
- Do not attempt to remove the breast shield from your breast while pumping. Turn the breast pump off and break the seal between your breast and breast shield with your finger, then remove breast shield from your breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- Contact your health care professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful. See Section 10 & 11 for more information.
- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain. For assistance with correct breast shield sizing and comfort please visit www.MedelaBreastshields.com or see a lactation consultant / breastfeeding specialist.
- Do not try to express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Make sure tubing is not kinked or pinched while pumping.
- Do not hold the pump kit by the bottle. This can lead to blockage of the milk ducts and engorgement.
- Using a breast pump on an aircraft is not recommended. The cabin pressure may affect the breast pumps performance.
- Wash hands thoroughly with soap and water before touching breast pump, kit and breasts, and avoid touching the inside of bottles or lids.
- Separate and wash all parts that are exposed to breast milk immediately after use. This will help remove breast milk residue and prevent growth of bacteria.

CAUTION: Can lead to minor injury: (cont.)

- Always inspect breast shields, connectors, valves, membranes, bottles, lids, and tubing prior to use for cleanliness. Contact Medela Customer Service if cleaning does not resolve the issue.
- Only use drinking-quality tap or bottled water for cleaning your breast pump and parts.
- Do not store wet or damp parts as mold may develop.
- Do not run pump with wet tubing. Doing so may damage the breast pump.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue against the breast shield tunnel, use of a lubricant such as Tender Care™ Lanolin may be beneficial. For assistance with correct breast shield sizing and comfort please visit www.MedelaBreastshields.com or see a lactation consultant/breastfeeding specialist.

SAVE THESE INSTRUCTIONS

Note

- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Bottles and component parts may become damaged if mishandled, e.g. dropped, over tightened, or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breast milk if bottles or components become damaged.
- Do not use Lithium batteries in the Pump in Style Advanced battery pack.
- Do not remove the breast pump from the soft case / carry bag.

Table of Contents

Important Safeguards	2-5
Product Description / Parts	8-9
1. Getting Started	10
2. Cleaning Overview	11
3. Wash - Before first use & after each use	12-13
3.1 Before using for the first time	12
3.2 After each use	13
4. Sanitize - Before first use & once per day	14-15
5. Breast Pump & Tubing Care	16-19
5.1 Tubing & faceplate care	16-19
5.2 Cleaning your motor unit	19
6. Assembly	20-27
6.1 Assembly of pump kit	20-23
6.2 PersonalFit™ Sizing	24
6.3 Powering your pump	25-27
6.3.1 Power adaptor care	25
6.3.2 Battery pack	26-27
6.3.3 Traveling outside of your country	26-27
7. Pumping	28-31
8. Storing Breast Milk	32
8.1 Storing breast milk in your Pump In Style® Advanced cooler bag	32
8.2 Storage Guidelines for Freshly Expressed Breast Milk	32
9. Preparing & Feeding Breast Milk	33
9.1 Preparing breast milk	33
9.2 Feeding breast milk.....	33
10.Troubleshooting	34
11.Supplemental Information	35
12.Disposal	36
13.Warranty	37
14.Meaning of Symbols	38
15.EMC/Technical Description	39-43
16.Technical Specifications	44

Take a moment to read through this entire instruction manual before using this product for the first time.

PLEASE SAVE THESE INSTRUCTIONS.

It is best to wait until your breastfeeding routine is established (approximately 4 weeks) before expressing breast milk, unless otherwise advised by your healthcare professional.

Glossary of Terms

2-Phase Expression® technology – research based technology that mimics a baby's natural nursing rhythm.

Expression Phase – slower sucking/pumping rhythm for gentle and efficient milk removal as quickly as possible.

Maximum Comfort Vacuum™ – highest vacuum setting where a mother still feels comfortable during pumping. It is different for every mother.

Stimulation Phase – fast sucking/pumping rhythm to stimulate the milk ejection reflex and to start the milk flowing.

Indications for use:

The Pump In Style Advanced breast pump is a powered breast pump to be used by lactating women to express and collect milk from their breast. The Pump In Style Advanced is a single user device.

Contraindications for use:

There are no known contraindications for use with this product.

i Tip

It is best to wait until your breast milk routine is established (approximately 4 weeks) before expressing breast milk, unless otherwise advised by your healthcare professional.

If you have medical reasons or other needs for exclusively pumping, it is recommended that you use a hospital-grade (multi-user) breast pump such as the Medela Symphony® Breast Pump. To learn more, visit www.medelabreastfeedingus.com.

Product Description

This breast pump is a personal-use electric breast pump that includes 2-Phase Expression® technology and is capable of single and double pumping.

The operating life of this breast pump is defined to be approximately three 15-minute sessions per day, for one year. The operating life for the breast pump kit is 6 months.

Product Description / Parts

Medela reserves the right to substitute any component or accessory with a replacement of equivalent performance.

2x PersonalFit™ Breast Shields, 24 mm

Applied Part (for other breast shield sizes see **Section 6.2**)

Replacement 87073

2x Connectors

Item 87071 or 8107264

87071


8107264

2x Valves

Item 87089

2x Membranes

Item 87088

4x Bottles

Item 6100050

Replacement 87130

4x Lids

Item 8107183

Replacement 87165

1x Battery pack

Item 9017002

Replacement 67553

2x Tubing

Item 87212

1x Power adaptor

Item 101036149

Replacement
101038300 or 68030

1x Cooler bag

Replacement 67068
(includes ice pack &
4 bottles with lids)

1x Ice pack

Item 8117010

Replacement 87092


On-the-go Tote

Backpack

1x Faceplate
Item 6007132

min
off
max

On/Off knob
Increase/reduce vacuum

pump in style
advanced

Let-down
button

Power
adaptor port

Tubing port Port plug Tubing port

UL US LISTED
SPS1 BREASTPUMP
9 VDC max.
Use Medela Item #9307910
Battery pack #907500

Additional parts included with The Metro Bag™

Pump In Style® Advanced
The Metro Bag™

27 mm PersonalFit™ Breast shields..... 87274

Metro Bag including:
Flexible, Plastic Workspace 3007299
2 Mesh Storage Bags

1. Getting Started

Breast shield assemblies come assembled. Separate all parts that will come into contact with the breast and breast milk before cleaning.

Parts to clean:

- Breast shields
- Breast milk bottles
- Lids
- Connectors
- Membranes
- Back caps

Note

It's important for you to do the following before using your Pump in Style Advanced Breast Pump for the first time:

- Separate all parts
- Clean - see instructions:
 - Wash - Section 3
 - Sanitize - Section 4

2. Cleaning Overview

When to Wash	Breast pump kit 	Breast pump bottles 
Before 1st use	<ul style="list-style-type: none"> ✓ wash ✓ sanitize 	<ul style="list-style-type: none"> ✓ wash ✓ sanitize
After each use	<ul style="list-style-type: none"> ✓ wash 	<ul style="list-style-type: none"> ✓ wash
Once per day	<ul style="list-style-type: none"> ✓ sanitize 	<ul style="list-style-type: none"> ✓ sanitize

When to Wash	Tubing 	Motor unit 
As needed	<ul style="list-style-type: none"> ✓ wash only if dirty or milk is present in tubing 	<ul style="list-style-type: none"> ✓ wipe with clean, damp cloth

Detailed cleaning instructions:

- Wash – Section 3
- Sanitize – Section 4
- Breast pump & tubing care – Section 5

Washing and sanitizing are two different activities and must be done separately to protect you and the performance of your breast pump.

3. Wash - Before first use & after each use

3.1 Before using for the first time

Clean and sanitize the parts below before using your pump for the first time. Follow the cleaning and sanitizing instructions in **Sections 3.2 and 4** that follow.

Supplies needed:

- Mild dish soap
- Clean dish cloth or soft brush
- Clean wash basin
- Drinking quality water

Parts to wash:

- Breast shields
- Breast shield connectors
- Membranes
- Valves
- Breast milk bottles and lids

Caution

- Only use drinking-quality tap or bottled water for cleaning.
- Take apart and wash all parts that come in contact with the breast and breast milk immediately after use to avoid dry up of milk residues and to prevent growth of bacteria.

3.2 After each use


Clean all parts that come in contact with your breast and breast milk. See image in **Section 1**.

Separate breast shields from breast shield connectors.

- Rinse all separated parts that come into contact with the breast and breast milk in cool water to remove breast milk residue.
- Once rinsed, parts that have touched the breast or breast milk should be placed in a clean wash basin, dedicated to only cleaning these items. Do not place parts directly in the sink.

Twist and pull the yellow valves off the breast shield connectors.

Remove the white membranes from the yellow valves.

Washing Options:

Washing in a Dedicated Basin

- Soak with warm soapy water for 5 minutes, wash with a clean, unused dish cloth and rinse with clean water.
- Place parts on a clean surface and/ or towel. Allow parts to air dry. Store dry parts in a cool place when not in use. Do not store wet or damp parts.

or

- You may wash your kit parts by following the instructions on the Quick Clean™ Breast Milk Removal Soap bottle. (Sold separately.)
- Wash all separated parts on top rack of dishwasher.
- Allow all breast pump parts to air dry in a clean area.
- Store dry parts when not in use.

4. Sanitize - Before first use & once per day

Note

- **Wash** breast pump kit parts after each use.
- **Sanitize** clean breast pump kit parts once daily.
- Take care not to damage parts of the pump kit during cleaning.
- If the individual components of the pump kit are cleaned in a dishwasher parts may be discolored by food pigments. This has no impact on the function.
- Store the pump kit in a clean bag/container until next use.

Supplies needed:

- Clean pot for boiling water
- Clean dish cloth
- Tongs
- Drinking quality water

Parts to sanitize:

- Breast shields
- Breast shield connectors
- Membranes
- Valves
- Breast milk bottles and lids

Sanitize daily

1

Refer to image in Section 1

Separate all clean parts that come in contact with your breast and breast milk.

4

Allow water to cool and gently remove parts from water with tongs.

Wash hands thoroughly.

- Let parts boil in water for 10 minutes.

or

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts when not in use. Do NOT store wet or damp parts.

You may sanitize your kit parts by following the instructions on Medela's Quick Clean™ Micro-Steam™ bag. (Sold separately.)

(Do not clean tubing in a micro-steam bag)

5. Breast Pump & Tubing Care

Caution

- Do not store wet or damp parts as mold may develop.
- If tubing becomes moldy, discontinue use and replace tubing. To find replacement parts visit www.ShopMedela.com.

Note

- Tubing should be washed if dirty or milk is present. Follow instructions in **Section 5.1**.
- Cleaning tubing is not necessary if condensation inside the tube is from previous washings or atmospheric conditions. Follow drying instructions in **Section 5.1**.

Supplies needed:

- Mild dish soap
- Clean dish cloth
- Clean wash basin
- Drinking quality water


Parts needed:

- Tubing
- Faceplate
- Diaphragm
- Breast pump bag


You may wash tubing by following the instructions on the Quick Clean™ Breast Milk Removal Soap bottle. (Sold separately.)

5.1 Tubing & faceplate care

Turn off breast pump.


Remove tubing from breast shield.


Info

Common causes of milk overflow into tubing:

- Heavy let-down
- Leaning forward or lying down while pumping
- Improper cleaning of valves and membranes
- Damaged membranes
- Overfilling bottles while pumping


2
Unplug breast pump from power source.

3
Remove tubing by pulling straight out of tubing ports. Do not wiggle or pull tubing at an angle.

5
Rinse tubing in cool water to remove breast milk.

6
Continued next page.

- Wash tubing in warm soapy water.
- Rinse tubing with clear water.
- Hang to air dry or connect the tubes to the running breast pump and allow any moisture that has formed to dry.


NOTE:

Clean the faceplate monthly. Also, inspect tubing after each pumping session. Clean tubing, faceplate and diaphragm if there is any sign of breast milk.

5. Breast Pump & Tubing Care (cont.)

Note

- Make sure faceplate and diaphragm are completely dry before reseating.
- The faceplate must be seated correctly and snapped together at all connection points. This ensures the pump will have adequate vacuum and optimum performance.

Snap off the faceplate using the thumb tab.

Locate the 3 connection points on the reverse side of the faceplate and on the pump.


- Soak in warm soapy water for 5 minutes.
- Clean with a clean dish-cloth or soft brush
- Rinse parts with clear water.
- Allow to air dry in a clean area.

- Some breast pump models may have a diaphragm that appears different than pictured above.
- Carefully wipe the entire surface area of the diaphragm where particles or liquid may have accumulated, with a clean (no soap), damp cloth.
- Allow all parts including edges to air dry, **completely**, in a clean area before replacing the faceplate.

5.2 Cleaning your motor unit

Align connection points and hold the back of the pump with one hand while pushing on the faceplate until it snaps into place.

Use a damp (not wet) cloth to wipe the breast pump bag after each use.

6. Assembly

Parts needed:

- Breast shields
- Breast shield connectors
- Membranes
- Valves
- Breast milk bottles and lids
- Tubing
- Breast pump
- 9 volt power adaptor

6.1 Assembly of pump kit

Wash hands before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.

Note


- Inspect valves for cracks or chips and inspect membranes for holes, tears or warping. If you notice any damage to these parts, replace them immediately.

Caution

- Wash hands thoroughly with soap and water before touching the breast and pump parts.


Info

- Carry out all the steps with care and assemble the pump kit correctly for optimum performance.

Snap membranes onto valves until membranes lie completely flat on valve.

2

Place your breast pump on a clean and dry surface. Ensure that the surface is solid and secure so your breast pump will not tip over or fall.

3

Securely attach breast shield to breast shield connector.

5

Push assembled yellow valves and white membranes onto the bottom of the breast shield connectors.

6

Screw bottles onto the breast shield connectors.

Continued next page.

6. Assembly (cont.)


Warning

- Use only the power adaptor supplied with the Pump In Style® Advanced.
- Always inspect the 9 volt power adaptor and Battery Pack cords prior to use for damaged, kinked, twisted, frayed or exposed wires. If damaged, immediately discontinue use. For replacement parts visit www.shopmedela.com.

Twists, kinks or bends

Has bulges

Exposed wires


Note

- See **Section 6.3.2** Battery Pack Assembly for instructions on the use of the Battery Pack.
- See **Section 6.3.3** Traveling Outside the U.S. for information on using your device internationally.

6.1 Assembly of pump kit (cont.)

7


87071


8107264

Insert the hard yellow adapter ends of tubing into the hole on the back of the breast shield connectors.

10

Plug the metal fitting on the end of the 9 volt power adaptor or Battery Pack cord into the power jack on the front of the breast pump faceplate.

8
Double pumping: Insert the free ends of tubing onto both ports on the front of the breast pump faceplate.

9
Single pumping: Insert the free end of tubing onto one port on the front of the breast pump faceplate. Place the port cap over the unused port.

11
Refer to Power Adaptor Care in Section 6.3.1

Plug the 9 volt power adaptor into a standard electrical outlet.

6. Assembly (cont.)

6.2 PersonalFit™ Sizing

1. Medela's PersonalFit sizing

Pumping should not hurt. For maximum comfort and pumping efficiency, Medela offers five breast shield sizes.

21
mm

24
mm

27
mm


30
mm

36
mm


Visit MedelaBreastshields.com to view the full guide to help you determine your optimal size based on your nipple diameter.

2. Test your breast shield size


1. Start with the 24 mm that came with your pump, or the size determined by measuring.
2. Center nipple and gently hold breast shield against your breast.
3. Adjust for Maximum Comfort Vacuum™, to achieve optimum suction level.
4. Refer to images while pumping in expression.

- Nipple rubs along side of tunnel
- Try a larger size

- Nipple is centered and moves freely

- Nipple and excessive areola are pulled into the tunnel
- Try a smaller size

3. Should you try a new size?

- Does your nipple rub sides of tunnel, to the point of causing discomfort?
- Do you see excessive areola being pulled into tunnel?
- Do you see any redness?
- Is your nipple or areola turning white?
- Do you feel unexpressed milk after pumping?

If you answered **“YES”** to any of these questions, consider trying a new size by following the measuring instructions above.

If you are still unsure if you selected the correct size, see a lactation consultant, breastfeeding specialist, or visit MedelaBreastshields.com for assistance in choosing the right size breast shield.

Caution

- While some discomfort may be felt when first using a breast pump, using the breast pump should not cause pain.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue with the breast shield tunnel, use of a lubricant such as Tender Care™ Lanolin may be beneficial.

6.3 Powering your pump

6.3.1 Power adaptor care


To care for your adaptor be sure to follow these storage directions.

Do not wrap the power adaptor cord around the plug body.

Correct!

Do not unplug the power adaptor by pulling on the cord.

Correct!

6. Assembly (cont.)

Note

- Remove the batteries if the Battery Pack is not going to be used for an extended period of time.
- If breast pump does not turn on, make sure the batteries have been inserted properly. If this does not solve the problem refer to the troubleshooting section of this instruction booklet or contact Medela Customer Service at 1-800-435-8316.
- Each set of batteries will provide approximately 2 hours of pumping time.
- Do not mix old and new batteries.
- Do not mix alkaline, standard (carbon-zinc) or rechargeable batteries.
- Do not use Lithium batteries.
- Only use alkaline or Ni-MH batteries
- For information about traveling with your breast pump and expressed breast milk, visit www.TSA.gov.

Replace battery cover.

6.3 Powering your pump (cont.)

6.3.2 Battery pack

You will need a total of 8 AA rechargeable batteries or alkaline batteries.

6.3.3 Traveling outside of your country

When traveling internationally, there are 3 power options:

or

Use the Battery Pack that was included with your Pump In Style® Advanced breast pump.

2

Replace all batteries at the same time.

Slide off both battery covers. Four (4) batteries are placed on each side of battery pack.

3

Place the **NEGATIVE (-)** end of the battery against the front of the spring by pushing the battery in and down in one motion until the battery snaps into place.

The **POSITIVE (+)** end of the battery should be near the positive sign on the battery pack.

or

Purchase the 9 volt Portable Vehicle Adaptor (sold separately) item 67174. Visit our online store at www.shopmedela.com.

The Portable Vehicle Adaptor may not work with all vehicles. Please check the owner's manual for compatibility with 12-24 volt input.

Use the power adaptor that was included with your Pump in Style® Advanced breast pump with universal prong converters (sold separately, not a Medela item)

Please note: The previous Pump in Style Advanced power adaptor (item number 68030 or 9207010) is not compatible outside of North America. Please call Medela Customer Service if you have any questions.

7. Pumping

Caution

- Always inspect breast shields, breast shield connectors, white membranes, yellow valves, and tubing prior to use for cleanliness. If dirty, see **Section 5.1**.
- If you experience pain or other symptoms while pumping seek advice from your healthcare professional or breastfeeding specialist. See **Section 11** for more information.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma may decrease milk output.
- Make sure tubing is not kinked or pinched while pumping.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk, or if expression is painful.

Info

- Research has shown that you will pump more efficiently – get more milk in less time – when pumping at Maximum Comfort Vacuum™ during the Expression Phase.
- You should reassess your Maximum Comfort Vacuum throughout your pumping experience, because it can change during the different stages of lactation, **see step 5 in Section 7**.
- Do not tilt or overfill containers when pumping.

Device operation

1

Wash your hands before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.

4

After two minutes, the pump will automatically change to the Expression Phase.

- a. If you let-down (your milk starts to flow) sooner than 2 minutes, press the let-down  button.
- b. If let-down has not happened during the Stimulation Phase and the pump has switched to the Expression Phase, press the let-down  button once to return to the Stimulation Phase.


2 While sitting in an upright position, center the assembled breast shields over your nipples.

3 To begin pumping, turn dial clockwise from the off position. Your pumping session will start in the Stimulation Phase. While pumping, adjust the speed/vacuum by turning the dial to your comfort level.

5 To find your Maximum Comfort Vacuum, increase speed/vacuum until pumping feels slightly uncomfortable (not painful), then decrease slightly.

6 When your pumping session is over, unplug tubing from the back of the breast shield(s) before setting down the bottles.

Continued next page.

7. Pumping (cont.)

Note

- Always disconnect Pump In Style® from the power source after expression.
- Only fill the bottle to the 150 mL mark.

Info

- **How often should you pump?** A breast pump is a replacement for when you are separated from your baby. It is important to pump when the baby would be breastfeeding, which usually means a working mother pumps 3 times during an 8 hour working day.
- **How long should your pumping session last?** Pumping times can vary from mother to mother sometimes 15 minutes, sometimes up to 20 minutes.
- **How should your breasts feel after pumping?** After pumping, your breasts should feel soft and there should be no firm areas. Before pumping, your breasts will have a firm, heavy feeling. After pumping make sure your breasts are not firm or have lumpy areas. This could indicate that the breast is not draining all over. See **Section 11** for more information.

Device operation (cont.)


7

Close the bottle with a lid.

- Follow the instructions in sections 8 and 9.

10

Unplug the breast pump from the wall outlet (power source)


8

Allow the pump to continue to run while you store your breast milk. This will allow any condensation that has formed from the natural results of humidity to dry in the tubing.

9

Turn off the breast pump.

11

Prior to storing your pump, remove tubing by pulling straight out of tubing ports. Do not wiggle or pull tubing at an angle.

12

Disassemble and clean following instructions in **Section 3**.

8. Storing Breast Milk

Warning

- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)

Note

- If the ice pack is frozen solid, it will safely cool your breast milk for up to 10 hours in a room temperature environment inside your pump bag or up to 8 hours in a room temperature environment, inside the cooler, stored outside of the pump bag.

8.1 Storing breast milk in your Pump In Style® Advanced cooler bag

The Pump In Style Advanced cooler bag requires one contoured ice pack (#87092, included) for cooling up to four 5 oz bottles of breast milk.

- As soon as you are finished pumping, place the breast milk bottles or breast milk storage bags in the cooler bag with the ice pack.
- Transfer breast milk bottles or breast milk storage bags to the refrigerator or freezer once you are home.

8.2

Storage Guidelines for Freshly Expressed Breast Milk (For Healthy Term Babies)

Room Temperature 60°F to 77°F (16° C to 25° C)	Refrigerator 39°F (4°C) or colder	Freezer 0°F (-18°C) or colder	Breast Milk Thawed in the Refrigerator
Up to 4 hours is best *Up to 6 hours for milk expressed under very clean conditions	Up to 3 days is best *Up to 5 days for milk expressed under very clean conditions	Up to 6 months is best *Up to 9 months for milk expressed under very clean conditions	At room temperature: Up to 2 hours Refrigerator: Up to 24 hours DO NOT REFREEZE

References: www.BreastmilkGuidelines.com

9. Preparing & Feeding Breast Milk

9.1 Preparing breast milk

- Thaw breast milk overnight in the refrigerator. Thawed breast milk is safe in the refrigerator for 24 hours. Do not refreeze thawed breast milk.
- Quickly thaw breast milk by holding the bottle under warm running water.
- Place the sealed bottle in a bowl of warm water for 20 minutes to bring it to body temperature. Make sure water does not come up to the level of the collar of the bottle.

9.2 Feeding breast milk

It is recommended that breastfeeding is well established prior to bottle feeding your baby.

- Always inspect the bottle, nipple and other components immediately before and after each use. If nipple appears cracked or torn, discontinue use immediately.
- To prevent possible choking hazard, test strength of nipple by pulling on bulb portion of the nipple.
- Do not attempt to enlarge the nipple hole.
- Infants must not be bottle fed without adult supervision.
- Nipple should not be used as a pacifier.

For additional breast milk collection and storage information, please visit:

www.BreastmilkGuidelines.com

10. Troubleshooting

<p>Low or No Suction</p>	<p>Valve and Membrane Check</p> <ul style="list-style-type: none"> • Remove and separate valve from membrane • Inspect the valve for chips or cracks. Inspect the membrane for holes or tears. If flawed or damaged replace valve and/or membrane prior to pumping • Wash valve and membrane per Cleaning section • Reassemble. Make sure the membrane lays flat on the valve, and if it does not, replace the valve and/or membrane prior to pumping <p>Breast shield and Breast shield Connector Check</p> <ul style="list-style-type: none"> • Inspect breast shield and breast shield connector to make sure there are no cracks and they securely attach to each other <p>Faceplate Check</p> <ul style="list-style-type: none"> • Wash faceplate and backplate per section 5. • Reassemble faceplate and backplate, and reattach tubing to faceplate port <p>Connection Point Check</p> <p>Inspect following connection points to ensure attachment is secure</p> <ul style="list-style-type: none"> • Tubing adapter to breast shield connector • Tubing end to faceplate port • Faceplate to backplate • Test vacuum. If suction is not improved, contact Medela Customer Service at 1-800-435-8316.
<p>Audible air coming from port cap/faceplate</p>	<p> CAUTION: Do not attempt to plug hole in cap or port with additional material.</p> <ul style="list-style-type: none"> • There is a built-in air leak for vacuum control, some audible air is normal.
<p>Power source</p>	<ul style="list-style-type: none"> • If pump is not turning on with the power adaptor try the battery pack. • If the pump is not turning on with the battery pack, be sure the batteries have been inserted correctly. • If pump is still not working with the power adaptor or the battery pack, contact Medela Customer Service at 1-800-435-8316.
<p>Faceplate port broken/cracked</p>	<p> CAUTION: If you are feeling a loss of suction with the broken faceplate, do not use the pump until the new faceplate arrives.</p> <ul style="list-style-type: none"> • Contact Medela Customer Service at 1-800-435-8316 or order online at www.shopmedela.com.
<p>Condensation or milk in tubing</p>	<ul style="list-style-type: none"> • Refer to "Breast pump & Tubing Care" section 5.
<p>No "let-down" or breast milk expression</p>	<ul style="list-style-type: none"> • Ensure your breast pump is assembled correctly and that there is suction. • Relax and take a 10-15 minute break if let down is not achieved after 2 consecutive pumping sessions. • Consult with your healthcare provider or breastfeeding specialist if expression does not occur.

11. Supplemental Information

Supplemental Information

Following are some common breastfeeding related conditions. **If you experience any of these symptoms, contact a healthcare professional or breastfeeding specialist.**

	Symptom	Potential Causes
Engorgement	Breasts are hard and uncomfortable, possibly with reddened areas.	Milk not draining properly from breast. Milk can drain back into the tissue causing areas to swell and become tender.
Blocked/Plugged Ducts	Area of breast will look reddened and may be tender to touch.	Milk not being drained from a specific duct. The area becomes 'clogged' up and milk is then prevented from flowing.
Mastitis	It often follows on from engorgement but can just occur suddenly. The initial symptoms are similar to flu – fatigue, headache and muscle aches, fever and localized breast tenderness. It can be in one or even both breasts and needs immediate treatment.	A bacterial infection in the breast tissue which is often accompanied by cracked nipples.
A Healthcare Professional or Breastfeeding Specialist can give you guidance in breastfeeding your baby and using a breast pump. Their guidance is valuable in terms of the lifetime of health benefits for both you and your baby. A listing of Lactation Professionals are available at www.medelalocator.com or call 1-800-TELL-YOU (1-800-835-5968).		

If you have not resolved the problem with your breast pump or you have further questions, please contact Medela Customer Service at 1-800-435-8316 or email customer.service@medela.com.

12. Disposal

Disposal

The unit is made of various metal and plastics. Before disposal, the device is to be rendered unusable and it must not be disposed of as unsorted municipal waste in accordance with local regulations. Use your local return and collection system for waste electrical and electronic equipment (incl. batteries). Improper disposal may have harmful effects on the environment and on public health.

Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Recycle your breast pump

Through the Medela Recycles program, Medela LLC provides customers the option to properly recycle their breast pump after they've finished their breastfeeding journey. For information on how to recycle your pump visit www.medelarecycles.com.

This program is currently only available to U.S. residents. We only accept the return of electric breast pumps. All other breast pump kit components (connectors, breast shields, tubing, bottles, valves, membranes, and cooling elements) can be disposed of in your waste management recycling container at your home. Cooler bags and other soft-good items are not recyclable.

13. Warranty

Caution

This is a single user product. Use by more than one person may present a health risk and performance questions and voids the warranty.

This product is warranted by Medela LLC to the original retail purchaser to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and accessories) from the date of purchase. Warranty can only be claimed in the country of purchase. In the event of a defect, Medela will repair or, at Medela's option, replace this product, without charge for such replacement, parts or labor. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO A DURATION OF 1 YEAR FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

Before you make a claim under this warranty, it may save time and expense to call Medela Customer Service (toll free) at 1-800-435-8316. You may also call this number for additional information concerning this warranty. If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, prepaid, together with your dated bill of sale or other proof of purchase and a brief statement of the problem to the following address:

Medela LLC – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS

Medela Canada
4160 Sladeview Crescent Unit #8
Mississauga, ON L5L 0A1
ATTENTION: RETURNS

Call first for authorization number. Returns not accepted without an authorization number.

14. Meaning of Symbols

Glossary of Symbols

Operator's Manual; Operating Instructions - To indicate that the operating instructions should be considered when operating the device.

This symbol indicates the manufacturer.

This symbol indicates do not dispose the device together with unsorted municipal waste (in accordance with local regulations).

This symbol indicates compliance with international requirements for protection from electric shock. (Type BF applied parts.)

IP22

This symbol indicates the protection against entrance of solid foreign objects and against harmful effects due to the entrance of water.

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates Class II.

This symbol indicates that the power adaptor is for indoor use only.

BREAST PUMP AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014), IEC 60601-1-6 (2010) + AMD 1 (2013), IEC 62366 (2007) + AMD 1 (2014) and IEC 60601-1-11 (2015)

This symbol indicates alternating current.

This symbol indicates direct current.

Caution Sign

UL Listed Part

Ten Pao Trademark

NOM Certified

15. EMC/Technical Description

The breast pump needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the instructions for use. Portable and mobile RF communications can affect the breast pump.

Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

15. EMC/Technical Description (cont.)

Guidance and manufacturer's declaration – electromagnetic emissions

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Emission tests	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Pst < 1.0	

Warning

Warning – This breast pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this breast pump should be observed to verify normal operation in the configuration in which it will be used.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 2kV, +/- 4kV, +/- 6 kV, +/- 8kV contact discharge +/- 2kV, +/-4kV, +/- 6 kV, +/-8kV, +/- 15 air discharge	+/- 2kV, +/-4kV, +/- 6 kV, +/- 8kV contact discharge +/- 2kV, +/-4kV, +/- 6 kV, +/-8kV, +/- 15 air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV , ± 1 kV Line-to-line	± 0,5 kV , ± 1 kV Line-to-line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	0 % UT ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the breast pump requires continued operation during power mains interruptions, it is recommended that the breast pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50 or 60 Hz	30 A/m, 50/60 Hz	It may be necessary to position the breast pump further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

15. EMC/Technical Description (cont.)

Guidance and manufacturer's declaration – electromagnetic immunity

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	6V	Portable and mobile RF communications equipment should be used no closer to any part of the breast pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommend separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz - 800 MHz $d = 2,3\sqrt{P}$ 800 MHz - 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b
Radiated RF IEC 61000-4-3	10 V/m (minimum) 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m (minimum)	Interference may occur in the vicinity of equipment marked with the following symbol. 

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the breast pump is used exceeds the applicable RF compliance level above, the breast pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the breast pump.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the breast pump

This electric breast pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the electric breast pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the electric breast pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2.3\sqrt{P}$
0,01	0,06	0,04	0,07
0,1	0,18	0,11	0,22
1	0,58	0,35	0,22
10	1,84	1,11	2,21
100	5,83	3,50	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the breast pump including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

16. Technical Specifications

Technical Specifications

Vacuum Range
-50.....-240 mmHg
54 ... 120 cpm

Operation
Temperature

Mode of operation: Continuous

Size
165 x 165 x 135 mm

Transport / Storage
Temperature

Weight
1.0 kg

Power adapter P/N - 101036149	
Power In	Power Out
100-240V~ 50/60 Hz 0.7A max	9.0 VDC 2A -----

Operation
Humidity

Battery Pack P/N - 9017002	
Power Out	
9.6 - 12 VDC 2A	

Transport /
Storage Humidity

Portable vehicle adaptor P/N - 9287011	
Power In	Power Out
12 -24 VDC 1.25A	9.0 VDC 1A -----

----- Direct current = "dc"

This symbol indicates compliance with international requirements for protection from electric shock. (Type BF applied parts.)

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Printed in the USA. Medela wordmark and logo, Invented by Medela 2-Phase Expression, Symphony, and Pump In Style are registered in the U.S. Patent and Trademark Office. Quick Clean, Micro-Steam, Maximum Comfort Vacuum, The Metro Bag and PersonalFit are trademarks of Medela.
Medela Canada, Inc.,